SPEAR TIP pH ELECTRODE

Model : PE-06HDA

ISO-9001, CE, IEC1010







LUTRON ELECTRONIC

The Art of Measurement

SPEAR TIP PH ELECTRODE

Model : PE-06HDA

FEATURES

The "Spear Tip pH electrode " is perfect for those pH measurements in

applications where sample piercing is required. Meat, sausage and cheese

are ideal applications. The electrode features a very durable glass

measuring spear packaged in a rugged virtually unbreakable epoxy body.

SPECIFICATIONS	
Measuring Range	1 to 13 pH (typical 0 to 14 pH).
Operation Temperature	5 to 60 C degree (41 to 140 F degree)
Electrode Structure	Combination type.
Reference of pH Electrode	Ag/AgCl, Double Junction.
Response time	\leq 2 minutes.
Body Material	Ероху.
Connector	BNC.
pH value of zero mV	7±1 pH.
Repetability	≧0.05 pH.
Mechanical Protection	With protection bottle on the electrode head.
Weight	68 g (including cable & socket).
Dimensions	Body length - 150 mm.
	Body Dia 12 mm.
	Cable length - 750 mm.

* Appearance and specifications listed in this brochure are subject to change without notice.

1612-PE06HDA